

**L-Spine Pre-Arthroplasty**

Siemens Flash

Application Examples: back pain

Oral Contrast	No
IV Contrast / Volume	No

*Technical Factors*

Detector Collimator	Acq 128 x 0.6 mm
Care kV	On / 120 kV
Care Dose 4D	On / 260 mAs
Rotation Time (seconds)	1.0
Pitch	0.8
Typical CTDIvol	17.55 mGy $\pm$ 50%

Topogram: AP and Lateral, 512 mm

Spine	Recon Type	Width / Increment	Algorithm	Safire	Window	FoV	Series Description	Networking	Post Processing
<b>Recon 1</b>	Axial	3 x 3	I70h	1	Bone	200	AXIAL BONE	PACS	None
<b>Recon 2</b>	Axial	3 x 3	I30s	1	Spine	200	AXIAL STND	PACS	None
<b>Recon 3</b>	3D:AXIAL	2 x 2	I70h	1	Bone	-	AXIAL MPR	PACS	Axial MPR
<b>Recon 4</b>	3D:COR	3 x 3	I30s	1	Bone	-	COR	PACS	Coronal MPR
<b>Recon 5</b>	3D:SAG	3 x 3	I30s	1	Bone	-	SAG	PACS	Sagittal MPR
<b>Recon 6</b>	Axial	0.6 x 0.6	I26s	1	Bone	200	AXIAL 0.6 STND	TeraRecon	None
<b>Recon 7</b>	3D:SAG	2 x 2	I30s	1	Bone	-	OBL RT	PACS	OBL MPR
<b>Recon 8</b>	3D:SAG	2 x 2	I30s	1	Bone	-	OBL LT	PACS	OBL MPR


This protocol is used for pre-operative planning for lumbar spine disc arthroplasty surgery. The images note the position of the pubic symphysis in relation to the disc space(s) to be replaced.

**Patient Position:** Supine feet first with arms comfortably above head.


**Scan Range:** One vertebral level above area of interest through pubic symphysis.

**Recons & Reformats:** Coronal FoV to include aorta through spinal anatomy and sagittal FoV to include the pubic symphysis anteriorly. Right and left oblique MPRs include facet joints.


1. Coronal MPR  
(through aorta)  
↓  
3x3mm




2. Sagittal MPR  
(large FoV)  
3x3mm




3. Axial MPR  
(through disc spaces)  
2x2mm



4. Oblique RT MPR  
(through right facet joints)  
2x2mm



5. Oblique LT MPR  
(through left facet joints)  
2x2mm



**3D:** No