

Child Craniosynostosis

Siemens 16 Slice

| |
|---|
| Application Examples: suture evaluation |
|---|

| | |
|----------------------|----|
| Oral Contrast | No |
| IV Contrast / Volume | No |

Technical Factors

| | |
|------------------------------------|-----------------------|
| Scan Type | Spiral |
| Detector Collimator | Acq 16 x 0.6 |
| kV / mAs / Rotation Time (seconds) | 110kV / 150 mAs / 1.0 |
| Care Dose 4D | Off |
| Pitch | 0.8 |
| Typical CTDIvol | 26.63 mGy |

Topogram: Lateral, 256 mm

| Head | Width / Increment | Kernel | Window | FOV | Series Description | Networking |
|---------|-------------------|--------|-------------|-----|-------------------------|------------|
| Recon 1 | 3 x 3 | C60s | Bone | 200 | AXIAL BONE | PACS |
| Recon 2 | 5 x 5 | C30s | Cerebrum | 200 | AXIAL STND | PACS |
| Recon 3 | 0.75 x 0.5 | C60s | Bone | 200 | AXIAL 0.75 x 0.5 BONE | MPR |
| Recon 4 | 0.75 x 0.5 | C20s | Base Orbita | 200 | AXIAL 0.75 x 0.5 SMOOTH | TERARECON |

This protocol is used for routine craniosynostosis studies.

Patient Position: Position head so the GML is perpendicular to the table in a symmetrical position (no rotation or tilt). Note no gantry angle with spiral acquisition. Axial images should be acquired parallel to a line drawn from the base of the skull to the glabella.

Scan Instructions: Position the head carefully to avoid compromising the airway. Baby's body may need to be elevated with a sponge or blanket to assure the head is in isocenter of the gantry. Use sponges to immobilize head in the head holder.

Scan requirements: Baby must remain motionless for entire scan. If sedation is needed, exam must be done in La Crosse CT.

Scan Range: Skull base through vertex. Scan in caudocranial direction.

2D Reformations: Post processing done in 3D card.

| Series: Craniosynostosis | Reformat Type | Width / Increment | Window | Series Description | Networking |
|--------------------------|---------------|-------------------|--------|--------------------|------------|
| Recon 3 | Coronal MPR | 3 x 3 | BONE | COR | PACS |
| Recon 3 | Sagittal MPR | 3 x 3 | BONE | SAG | PACS |

3D: VR and skull views. Contact La Crosse Imaging Lab.