

Child Face

Siemens 16 Slice

Application Examples: facial trauma, fracture

Oral Contrast	No
IV Contrast / Volume	No

Technical Factors

Scan Type	Spiral
Detector Collimator	Acq 16 x 0.6 mm
kV / mAs / Rotation Time (seconds)	110 kV / 50 mAs / 1.0
Care Dose 4D	On
Pitch	0.8
Typical CTDIvol	8.88 mGy

Topogram: Lateral, 256 mm

Face	Width / Increment	Kernel	Window	FoV	Series Description	Networking
Recon 1	1 x 1	C60s	Bone	150	AXIAL BONE	PACS
Recon 2	1 x 1	C30s	Larynx	150	AXIAL STND	PACS
Recon 3	0.75 x 0.5	C60s	Bone	150	AXIAL 0.75 x 0.5 BONE	MPR / Definition
Recon 4	0.75 x 0.5	C20s	Bone	150	AXIAL 0.75 x 0.5 SMOOTH	TeraRecon

This protocol is used for routine facial bone studies.


Patient Position: Position patient so IOML is perpendicular to table and head is in a symmetrical position (no rotation or tilt).

Scan Range: Frontal sinus through hard palate. Include mandible if requested or if the injury is to the mandible.

Reformations: Post processing done in 3D card. Coronal and sagittal MPRs are reconstructed perpendicular to hard palate. Extend coronal MPR through cervical spine. If unable to place patient in ideal position, make axial MPR data set parallel to hard palate – 1mm x 1mm.

Series: Face	Reformat Type	Width / Increment	Window	Series Description	Networking
Recon 3	Coronal MPR	3 x 3	BONE	COR	PACS
Recon 3	Sagittal MPR	3 x 3	BONE	SAG	PACS


Coronal MPR
3x3mm



Sagittal MPR
3x3mm



Axial MPR
1x1mm



3D: Upon request.