

Shoulder

Siemens Flash

| |
|---|
| Application Examples: fracture, dislocation |
|---|

| | |
|----------------------|----|
| Oral Contrast | No |
| IV Contrast / Volume | No |

| | |
|-------------|-------------|
| Breath Hold | Inspiration |
|-------------|-------------|

Technical Factors

| | |
|-------------------------|---------------------|
| Detector Collimator | Acq 128 x 0.6 mm |
| Care kV | On / 120 kV |
| Care Dose 4D | On / 200 mAs |
| Rotation Time (seconds) | 1.0 |
| Pitch | 0.6 |
| Typical CTDIvol | 13.52 mGy \pm 50% |


Topogram: Lateral and AP, 256 mm

| Shoulder | Recon Type | Width / Increment | Algorithm | Safire | Window | FoV | Series Description | Networking | Post Processing |
|----------------|------------|-------------------|-----------|--------|----------|-----|--------------------|------------|-------------------|
| Recon 1 | Axial | 3 x 3 | I70h | 2 | Shoulder | 200 | AXIAL | PACS | None |
| Recon 2 | 3D:COR | 2 x 2 | I70h | 2 | Shoulder | - | COR | PACS | Coronal MPR |
| Recon 3 | 3D:SAG | 2 x 2 | I70h | 2 | Shoulder | - | SAG | PACS | Sagittal MPR |
| Recon 4 | 3D:AXIAL | 2 x 2 | I70h | 2 | Shoulder | 200 | OBL AXIAL | PACS | Oblique Axial MPR |
| Recon 5 | Axial | 0.6 x 0.6 | I26s | 2 | Shoulder | 200 | AXIAL 0.6 STND | TeraRecon | None |

Patient Position: Patient lying in supine position, head first, shoulders square with affected shoulder slightly toward iso-center. Affected arms should be in neutral rotation. Unaffected arm positioned above head on large patients.

Scan Range: Scan entire gleno-humeral joint and through area of interest. If for scapula, include entire scapula in scan range.

Recons and Reformations: Coronal, sagittal and oblique axial MPRs should be made in orthogonal planes to gleno-humeral joint as depicted below.



3D: Upon request. See post processing protocol.