

Abdomen Three Phase KUB

Siemens go.All

Application Examples: hematuria

Oral Contrast	H2O
IV Contrast / Volume <300lbs	150 ml Omnipaque 300
IV Contrast / Volume >300lbs	150 ml Omnipaque 350
Injection Rate	3ml / sec (no split bolus)

Technical Factors

Renal Calc & Delayed KUB	
Detector Collimator	Acq 32 X 0.7 mm
Care kV	On / 120 kV
Care Dose 4D	On / 110 mAs
Rotation Time (seconds)	0.5
Pitch	0.8
Typical CTDIvol	10.34 mGy

Liver/Kidneys	
Detector Collimator	Acq 32 X 0.7 mm
Care kV	On / 120 Kv
Care Dose 4D	On / 110 mAs
Rotation Time (seconds)	0.5
Pitch	0.8
Typical CTDIvol	10.34 mGy

Topogram: Lateral & AP, 512 mm

Renal Calc	Recon Type	Width / Increment	Algorithm	Safire	Window	Series Description	Networking	Post Processing
Recon 1	Axial	3 x 3	Br40	2	Abdomen	AXIAL WITHOUT	PACS	None

AP	Recon Type	Width Increment	Algorithm	Safire	Window	Series Description	Networking	Post Processing
Recon 1	Axial	3 x 3	Br40	2	Abdomen	AXIAL	PACS	None

Delayed	Recon Type	Width / Increment	Algorithm	Safire	Window	Series Description	Networking	Post Processing
Recon 1	Axial	3 x 3	Br40	2	Abdomen	AXIAL DELAYED	PACS	None
Recon 2	3D:COR	3 x 3	Br40	2	Abdomen	COR	PACS	Coronal MPR
Recon 3	3D:SAG	3 x 3	Br40	2	Abdomen	SAG	PACS	Sagittal MPR
Recon 4	3D:COR	3 x 3	Br40	2	Abdomen	COR MIP	PACS	Coronal MIP
Recon 5	Axial	0.6 x 0.6	Br36	2	Abdomen	AXIAL DELAYED 0.6 STND	TeraRecon	None

This protocol is used for evaluating common causes of persistent hematuria such as stones or tumors.

Exam Instructions: Patient should be instructed to drink one quart of water prior to arrival. If patient arrives without drinking water prior, give patient one quart of water to drink approximately 30 minutes before scan.

Patient Position: Patient lying supine with arms above head.

Scan Instructions: First, scan non-contrast kidneys through bladder. Second, inject IV contrast all 150ml and scan entire liver through bladder at 90 second scan delay. Third, wait 10 minutes and scan kidneys through bladder.

Recons and Reformations: Adjust FoV to fit body contour.

3D: Raysum

**IMAR on hematuria protocols should be used on the non contrast axials and delayed axial phase images. (only need through the area of metal)