

Chest without Contrast

Siemens go.All

Application Examples: Pulmonary nodule

Technical Factors


Breath Hold	Inspiration
Scan Type	Spiral
Detector Collimator	Acq 32 x 0.7 mm
X-Care	On
Care kV	On / 120 kV
Care Dose 4D	On / 60 mAs
Rotation Time (seconds)	0.5
Pitch	0.6
Typical CTDIvol	5.5 mGy ± 50%

Topogram: Lateral and AP, 512 mm

Chest	Recon Type	Width / Increment	Algorithm	Safire	Window	Series Description	Networking	Post Processing
Recon 1	Axial	5 x 5	Br40	2	Mediastinum	AXIAL	PACS	None
Recon 2	Axial	2 x 2	Br64	2	Lung	AXIAL LUNG	PACS	None
Recon 3	3D:COR	3 x 3	Br40	2	Mediastinum	COR	PACS	Coronal MPR
Recon 4	3D:SAG	3 x 3	Br40	2	Mediastinum	SAG	PACS	Sagittal MPR
Recon 5	3D: COR	8 x 5	Br40	2	Lung	COR MIP	PACS	Coronal MIP
Recon 6	3D: AXIAL	8 x 5	Br40	2	Lung	AXIAL MIP	PACS	Axial MIP
Recon 7	Axial	1.0 x 0.8	Br36	2	Mediastinum	AXIAL 1.0 x 0.8 STND	TeraRecon	None
Recon 8	Lung CAD	1.0 x 0.7	Br60	2	Lung	LUNG CAD	PACS	None

If pt is over 250lbs use CHEST WO over 250lbs on scanner.**Patient Position:** Patient lying supine with arms above head and lower legs supported.**Scan Range:** Lung apices through adrenal glands.**Recons and Reformations:** Set recon 2 begin and end points to include *lungs only*. If indication is for follow up aortic aneurysm, create oblique “candy cane” view.**Candy Cane View**—Copy the AXIAL recon box, change Image Orientation to Sag, Image order Rt to Lt, Image Type-MPR

Chest	Recon Type	Width / Increment	Algorithm	Safire	Window	Series Description	Networking	Post Processing
Recon 9	3D:OBL	3 x 3	Bv40	2	Mediastinum	OBL MPR	PACS	None



Set up

End Result

If being done for TRAUMA use the appropriate protocol labeled TRAUMA or add 2 extra Recon boxes

- 1) Radial Rib
- 2) Rib Parallel