

**Femoral Anteversion**

Siemens go.All

Application Examples: femoral anteversion

*Technical Factors*

| <b>Hips</b>              |                        |
|--------------------------|------------------------|
| Detector Collimator      | Acq 32 x 0.7mm         |
| kV / mAs / Rotation Time | 100 kV / 120 mAs       |
| Care Dose 4D             | Off                    |
| Pitch                    | 0.9                    |
| Typical CTDIvol          | 4.94 mGy ± 50%         |
| <b>Knees</b>             |                        |
| Detector Collimator      | Acq 32 x 0.7mm         |
| kV / mAs / Rotation Time | 100 kV / 100 mAs / 1.0 |
| Care Dose 4D             | Off                    |
| Pitch                    | 0.9                    |
| Typical CTDIvol          | 4.12 mGy ± 50%         |
| <b>Ankles</b>            |                        |
| Detector Collimator      | Acq 32 x 0.7mm         |
| kV / mAs / Rotation Time | 100 kV / 80mAs / 1.0   |
| Care Dose 4D             | Off                    |
| Pitch                    | 0.9                    |
| Typical CTDIvol          | 3.30 mGy ± 50%         |

Topogram: Lateral &amp; AP, 1024 mm


| Hips    | Recon Type | Width / Increment | Algorithm | Window | Series Description | Networking | Post Processing |
|---------|------------|-------------------|-----------|--------|--------------------|------------|-----------------|
| Recon 1 | Axial      | 3 x 3             | Br64      | BONE   | HIPS AXIAL         | PACS       | None            |

| Knees   | Recon Type | Width / Increment | Algorithm | Window | Series Description | Networking | Post Processing |
|---------|------------|-------------------|-----------|--------|--------------------|------------|-----------------|
| Recon 1 | Axial      | 3 x 3             | Br64      | BONE   | KNEES AXIAL        | PACS       | None            |

| Ankles  | Recon Type | Width / Increment | Algorithm | Window | Series Description | Networking | Post Processing |
|---------|------------|-------------------|-----------|--------|--------------------|------------|-----------------|
| Recon 1 | Axial      | 3 x 3             | Br64      | BONE   | ANKLES AXIAL       | PACS       | None            |

*Femoral Anteversion* is a condition in which the femoral neck leans forward with respect to the rest of the femur. This causes the lower extremity on the affected side to rotate internally (i.e. the knee and foot twists toward the midline of the body). The primary purpose of this scan is to allow the radiologist to measure the angle of rotation of the femoral necks relative to the femoral condyles bilaterally. Similar measurements can also be made of the tibias.

**Patient Position:** Supine, feet first, legs flat on the table (no cushions or wedges under legs or feet). Position legs as close together as possible in their **neutral position**. Consider taping feet and/or knees together to help stabilize if needed, but do not internally rotate legs.



**Topogram:** Lateral & AP views to include acetabulum through malleoli.

**Scan Instructions:** Acquire three separate spiral data sets keeping all FoVs and X and Y data points consistent. First set hips range (top of femoral head through lesser trochanters). Second set knees range (femoral condyles through tibial plateau or growth plate to growth plate). And third set to include ankles (growth plate through malleoli).

**Scan Note:** If patient is an adult with a prosthetic knee, scan coverage for knee is to include entire prosthesis.

**3D:** Yes. Notify 3D Technologist exam is ready for post processing.