

## Radio-embolization Liver

Siemens go.All

Application Examples: pre or post radio-embolization of liver	
Scan Protocol (See for Technical Factors)	<b>Liver</b>
Oral Contrast	<b>No</b>
IV Contrast / Volume	Omnipaque 350 / P3T

This protocol is used to evaluate pre or post radio-embolization of liver.

**Patient Position:** Patient lying in supine position, arms positioned comfortably above head and lower legs supported.

**Scan Instructions & Range:** Use 4 phase liver protocol. Image all phases.