

ELBOW-STRESS VIEW

ROUTINES: Elbow Stress View.

DISTANCE: 40"

BUCKY: No

FOCAL SPOT: Small

IMAGING PLATE: 8x10 CR cassette or DR detector

CENTRAL RAY: Perpendicular to elbow joint.

HELPFUL POSITIONING:

1. Patient lying on table with arm abducted away from body.
2. Hand externally rotated to bring elbow to AP position.
3. Use small table or cart to hold cassette in place.

Reviewed 10/01/19

Neutral x-table elbow

Stress view x-table

