

KNEE-TUNNEL (HOLMBLAD)

ROUTINE: Holmblad

DISTANCE: 40"

BUCKY: Yes

FOCAL SPOT: Small

IMAGING PLATE: 10x12 CR cassette or DR detector

CENTRAL RAY:

1. Perpendicular to film entering the popliteal region and exiting $\frac{1}{2}$ " inferior to the patellar apex.

HELPFUL POSITIONING:

1. Patient is kneeling on the table. Patient leans forward so that the affected knee is flexed 70° to the tabletop.
2. If patient cannot kneel take a *Camp Coventry* view instead (see pg 20).



Reviewed 10/01/19