

## PELVIS-JUDET VIEWS

**ROUTINE:** 45° Obliques of the Pelvis.

**DISTANCE:** 44"

**BUCKY:** Yes

**FOCAL SPOT:** Large

**IMAGING PLATE:** 14x17 CR cassette or DR detector

**CENTRAL RAY:**

Mid pelvis – perpendicular to film.

**HELPFUL POSITIONING:**

Patient supine, oblique patient 45° both ways (RPO & LPO).

**ADDITIONAL COMMENTS:**

Judet views of a single hip may be ordered if doctor desires.

*Reviewed 10/01/19*

