

## **SACRUM/COCCYX**

**ROUTINE:** AP of both the Sacrum & the Coccyx separately & Lateral of both on one.

**DISTANCE:** 40"

**BUCKY:** YES

**FOCAL SPOT:** Large

**IMAGING PLATE:** 10x12 CR cassettes or DR detector

### **CENTRAL RAY:**

***Sacrum:*** 15° cephalic angle entering at the midline halfway between the pubic symphysis and the ASIS.

***Coccyx:*** 10° caudal angle entering 2" superior to the pubic symphysis.

*Reviewed 10/01/19*