

MRI Enterography 3T Only

Application Examples: small bowel disease or crohns disease

Patient Prep: Small bowl prep day before. Brezza 900 mL, (one bottle (450 mL) every 20 minutes while in MR

Fov Coverage: Top of colonic flexures to anus

Sequence	FoV	Thickness	Gap	Gap
Coronal Haste/SSFSE	40.0 cm	5.0 mm	1.0 mm	Abdomen and Pelvis (AP)
Axial Haste/SSFSE	40.0 cm	6.0 mm	0.0 mm	AP
Axial Haste FS/SSFES FS	40.0 cm	6.0 mm	0.0 mm	AP
Inject Glucagon 1mL				
Coronal Vibe FS/Lava Flex pre	40.0	4.4mm	0.0	AP
INJECT CONTRAST (1mmol/Kg)				
Coronal Vibe FS/Lava Flex Post	40.0 cm	4.4 mm	0.0 mm	45 Sec after injection, AP
Coronal Vibe FS/Lava Flex Post	40.0 cm	4.4 mm	0.0 mm	2min after injection, AP
Axial Vibe FS/ Lava Flex Post	40.0 cm	4.4 mm	0.0 mm	AP