
Gundersen Health System Neuro MR Protocol GE 1.5T 450W

Thoracic Spine Application Examples: DJD, dorsal back pain

Sequence	FoV	Thickness	Gap
SAG T1	35CM	3MM	0.5MM
SAG T2	35CM	3MM	0.5MM
SAG STIR	35CM	ЗММ	0.5MM
AX T2 UPPER			
(MAKE SURE TO COVER ENTIRE TSPINE)	20CM	4MM	0.5MM
AX T2 LOWER			
(MAKE SURE TO COVER ENTIRE TSPINE)	20CM	4MM	0.5MM



SAG Coverage

AXIAL Coverage