

Upper Extremity Neurogram 3T ONLY

Application Examples: ulnar nerve numbness or pain

Split axials into two separate segments to cover area. May limit area if ordering provider specifies area to be scanned (for example: elbow or cubital tunnel).

Sequence	FoV	Thickness	Gap
AX T1	14CM	3MM	1MM
AX T2 IDEAL	14CM	3MM	1MM
COR STIR	15CM	3MM	1MM
INJECT CONTRAST (1mmol/kg)			
AX T1 FS WITH GAD	14CM	3MM	1MM
COR T1 FS WITH GAD	15CM	3MM	1MM