

## ***Lower Extremity Neurogram 3T ONLY***

Application Examples: nerve numbness or pain

Split axials into two separate segments to cover area. May limit area if ordering provider specifies area to be scanned (for example: knee or perineal nerve).

<b>Sequence</b>	<b>FoV</b>	<b>Thickness</b>	<b>Gap</b>
AX T1	14CM	3MM	1MM
AX T2 IDEAL	14CM	3MM	1MM
COR STIR	15CM	3MM	1MM
INJECT CONTRAST (1mmol/Kg)			
AX T1 FS WITH GAD	14CM	3MM	1MM
COR T1 FS WITH GAD	15CM	3MM	1MM