

Routine Ankle Arthrogram

Application Examples: Loose Body

Sequence	FoV	Thickness	Gap
Axial T1 FS	14.0	3.0mm	0.8mm
Axial T2 FS	14.0	3.0mm	0.8mm
Coronal T1 FS	14.0	3.0mm	0.8mm
Coronal T2 FS	14.0	3.0mm	0.8mm
Coronal T1 Non FS	14.0	3.0mm	0.8mm
Sagittal T2 FS	14.0	3.0mm	0.8mm
Sagittal T1 Non FS	14.0	3.0mm	0.8mm

