

Routine SI Joint

Application Examples: SI joint arthritis, SI joint disfunction

Sequence	FoV	Thickness	Gap
Sagittal T2 FS (L5-Coccyx)	28.0	4.0mm	1.0mm
Coronal Obl T1	20.0	3.0mm	1.0mm
Coronal Obl T1 FS	20.0	3.0mm	1.0mm
Coronal Obl Stir	20.0	3.0mm	1.0mm
Axial Obl T1	20.0	3.0mm	1.0mm
Axial Obl Stir	20.0	3.0mm	1.0mm
Inject Contrast (1.0mmol/KG)			
Axial Obl T1 FS	20.0	3.0mm	1.0mm
Coronal Obl T1 FS	20.0	3.0mm	1.0mm



Sagittal (L5-Coccyx) Coverage



Coronal Obl Coverage



Axial Obl Coverage