

Routine Finger

Application Examples: pain, fracture

Scan Range: First sequence (coronal PD) is always entire hand. For thumb imaging, align slices as depicted below and for fingers align straight to joint.

Sequence	FoV	Thickness	Gap
Coronal PD	9.0	2.0mm	0.5mm
Axial T1	9.0	2.0mm	0.5mm
Axial T2 FS	9.0	2.0mm	0.5mm
Coronal T1	9.0	2.0mm	0.5mm
Coronal T2 FS	9.0	2.0mm	0.5mm
Sagittal T1	9.0	2.0mm	0.5mm
Sagittal T2 FS	9.0	2.0mm	0.5mm
Sagittal Fiesta	10.0	2.0mm	0.5mm



Axial Coverage



Coronal Coverage



Sagittal Coverage