

IMG1200 (Adult)
IMG1265 (Child under 12)

PELVIS

ROUTINE: AP Pelvis.

DISTANCE: 44"

BUCKY: Yes

FOCAL SPOT: Large

IMAGING PLATE: 14x17 CR cassette or DR detector

CENTRAL RAY:

1. Perpendicular entering the midline 2" superior to the pubic symphysis.

HELPFUL POSITIONING:

1. Patient is supine with legs extended & internally rotated 15°.

ADDITIONAL COMMENTS:

1. **Pelvis for Hips**-center lower-top of film at ASIS. Be sure to include all the prosthesis. IMG1200
2. **Pelvis for Sacroiliitis**-do a PA pelvis. IMG1200

SEE ALSO/OTHER VIEWS:

1. Ferguson view of Pelvis-30° cephalic tube angle done to demonstrate L4-L5 & SI joints.

IMG1200 PELVIS FOR PAIN OR TRAUMA
IMG1205 PELVIS AND EITHER HIP
IMG1162L 2V LEFT HIP
IMG1162R 2V RIGHT HIP
IMG1107 PELVIS & BOTH LATERAL HIPS
IMG1265 PELVIS FOR INFANT & CHILD
IMG1265 PELVIS & FROGLEG PELVIS FOR INFANT & CHILD

Updated 03/03/2021 rsj